

Light Cure Bracket Adhesives



NEWLY DEVELOPED - SUPERIOR WORKING PROPERTIES

- Proven bond strength and consistency
- Unique viscosity easily penetrates bracket pads
- No bracket drift – stays in place until cured
- New proprietary microfiller construction
- One-piece flash removal – "ropes" off bracket periphery
- Translucent finish matches any tooth shade
- Maximum bond strength
- Easy removal and cleanup
- UV disclosing agent – fluoresces on the tooth under black light

CLINICAL APPLICATIONS

- Bracket and tube bonding
- Fixed lingual retainer bonding
- Bonds all metal, ceramic and composite attachments
- Reliamold compatible (see Reliamolds on page 29)
- Clear aligner attachment material
 - Viscous to allow for preloading of trays

GO TO®

4g Paste Syringe	#GTP
4g Paste In Tips (10 Tips)	#GTT



EXCESS ADHESIVE ROPES OFF IN A SINGLE PIECE FROM THE PERIPHERY OF THE BRACKET



BLACKLIGHT FLUORESCENT



CRYSTAL CLEAR FOR UNIVERSAL SHADE MATCHING WITH CLEAR ALIGNERS

RECOMMENDED PRIMER WITH
GOTO® ADHESIVE

Assure® Plus (Page 3)



- * Etchants and Conditioners appear on page 17
- * Bonding Accessories appear on page 35
- * Auxiliary Bonding Supplies appear on page 35



GO TO ADHESIVE WORKING PROPERTIES ARE IDEAL FOR CLEAR ALIGNER ATTACHMENTS



1. PUMICE



2A. ETCH



2B. ETCH



3. POST ETCH CHALKY SURFACE



4A. ASSURE PLUS



4B. ASSURE PLUS APPLICATION



4C. ASSURE PLUS AIR DRY



4D. ASSURE PLUS CURE
(OPTIONAL BUT RECOMMENDED)



5. ADHESIVE APPLICATION



6. TRAY SEATING



7. CURE



8. FUNCTIONAL ATTACHMENT