JUST SCAN IT

Aoralscan 3 An Intra Oral Scanner.





Technical Specifications

Aoralscan 3


Scan Field	Standard tip: 16 x 12 x 22 mm Small tip: 12 x 9 x 22 mm
Scan Principle	Non-contact with structured light
Dimension (LxWxH)	281 x 33 x 46 mm
Weight	240 ± 10 g (Without Cables)
Output	STL, OBJ, PLY
Connector	USB 3.0
Power Input	12V DC/3A

Recommended PC Configuration

Aoralscan 3

CPU	Intel Core i7-8700 or higher
RAM	16 Gigabytes (GB) or more
Hard Drive	256 Gigabytes (GB) SSD or more
GPU (Graphic Card)	NVIDIA® RTX 2060 6GB or better
Operating System	Windows 10 Prodessional (64-bit) or later versions of Windows operating system
Screen Resolution	1920 x 1080, 60Hz or more
I/O Ports	More than 2 type-A USB 3.0 (Or higher) ports







Capture Digital Impressions in Real Time

Fast, Accurate and Exceptional Scan Performance

The constant optimization and innovation of SHINING 3D's hardware and software system, led to the release of its latest and finest intraoral scanner, Aoralscan 3.



22 mm




scan bodies and periodontal scans.

Aoralscan 3 is suitable for various clinical scenarios including acquiring

Clinical Toolkit

The clinical toolkit assists dentists to evaluate and pre-design the scanned data in their own clinic. It facilitates effective communication with dental labs to deliver quality dental work with high efficiency.



User-Friendly

Smart and Intuitive - Smooth Scanning Experience



AI SCAN

Automatically identifies and filters out unnecessary soft tissue data during scanning resulting in a quicker and cleaner scan process.

MOTION SENSING

Dentists can scan with minimal contact with the computer, ensuring a safe and efficient scan experience.

Distributed By

